Stereospecific Detection Technologies Reagents beyond the chemistry

SDT GmbH • Pascalstraße 17 52499 Baesweiler • Germany www.sdt-reagents.de

Product Specification

Product Name: 5-Bromo-5-nitro-1,3-dioxane (BND), 0,5% aqueous solution

Product Code: #BND-W

Lot Number: 693164

Storage: room temperature (+18°C~+28°C). NOTE that that at temperatures below 18°C BND may crystallize out of the saturated solution. DO NOT STORE this product at temperatures lower than 18°C!

Expiration date: 07/2021

Appearance: Transparent non viscous liquid

Odour: Weak, characteristic

pH: 5,5 - 6,5 (with NO buffer capacity)

Preservative/anti-microbial: 5-Bromo-5-nitro-1,3-dioxane, 5000 ppm

<u>Description:</u> #BND-W is essentially a saturated 0,5% solution of 5-Bromo-5-nitro-1,3-dioxane (BND) in the pure deionized water, 0.2 µm filtered. BND is known as a well proven and widely validated low-hazardous IVD preservative giving dependable anti-microbial protection to the diverse stable liquid IVD product formulations.

Recommended use as example of making n ml buffer preserved with e.g. 0,1% BND:

- 1. Pour into a glass 1/5 n (gravimetrically measured as 1 ml ~ 1 g) of #BND-W.
- 2. Add 4/5 n of deionized water.
- 3. Add salts and/or other buffering reagents as per given buffer formula, mix, measure and adjust pH.
- 4. Filter through 0,2 µm filter unit.

This product is for *in vitro* research or further IVD manufacturing use only.